


How to use


1 Turn on



2 Select power level,
wait for the signal



3 Place your finger



Advantages of the Product

NO CONTACT

No infection

Minimal maintenance

Reduce pain


Low electrical consumption

Portable device

NO WASTE

What's new?

Onteko has developed a new lancet for diabetic patients. Using state of the art laser technology, it is the safest choice in the market.



What is it ?

The Laser Lancet has different energy levels, according to different types of skin. It detects the finger once it is on the focal point, then emits a painless laser pulse that perforates the skin. It doesn't have any contact with the skin, and, since it does not require consumable parts, there is no waste.



onteko

laser solutions

Contact Information

email

info@ontekollc.com

website

www.onteko.net



onteko

laser solutions

